

## ***Microsatellite Libraries by Enrichment.***

P A Bloor, F S Barker, P C Watts, H A Noyes, S J Kemp

Version 1.0, September 2001.

**Animal Genomics Laboratory,**  
Lab 1.03g,  
Donnan Laboratories,  
Crown Street,  
School of Biological Sciences,  
University of Liverpool,  
L69 7ZD. UK.

<http://www.liv.ac.uk/~kempsj/genomics.html>

This protocol was developed by Faye Barker and Paul Bloor as part of their respective PhD projects. This document was assembled by Phill Watts who has used it to guide a succession of undergraduate students with spectacular success. It derives from a mix of protocols described in the references on the next page. We use Qiagen kits and Promega vector cloning kits - we *strongly* recommend that you read these manuals before using the kits. Alternative kits may work just as well.

***We hope it works as well for you - If you publish work using this protocol, please cite it.***

Any comments or suggestions should be sent to Steve Kemp [ [kempsj@liv.ac.uk](mailto:kempsj@liv.ac.uk) ]

**References:**

- Fischer, D. & Bachmann, K. 1998. Microsatellite enrichment in organisms with large genomes (*Allium cepa* L.. *Biotechniques*. **24**, 796-802.
- Gardner, M.G., Cooper, S.J.B., Bull, C.M. & Grant, W.N. 1999. Isolation of microsatellite loci from a social lizard, *Egernia stokesii*, using a modified enrichment procedure. *Journal of Heredity*. **90**, 301-304.
- Hamilton, M.B., Pincus, E.L., Di Fiore, A. & Fleischer, R.C. 1999. Universal linker and ligation procedures for construction of genomic DNA libraries enriched for microsatellites. *Biotechniques*. **27**, 500-507.
- Refseth, U.H., Fangan, B.M. & Jakobsen, 1997. Hybridization capture of microsatellites directly from genomic DNA. *Electrophoresis*. **18**, 1519-1523.
- Sambrook, J., Fritsch, E.F. & Maniatis, T. 1989. Molecular cloning: a laboratory manual, 2nd ed. Cold Spring Harbour, New York: Cold Spring Harbour Laboratory Press.

**Materials required:**

Glucose,  
isopropanol, LA, LB,  
MgCl<sub>2</sub>.6H<sub>2</sub>O, MgSO<sub>4</sub>.7H<sub>2</sub>O,  
SOC-medium, 20 X SSC, 6 X SSC, 2 X SSC, 1 X SSC, Tris-HCl (10 mM, pH 7.5),  
1 X Washing/Binding (W/B) buffer (1 M NaCl, 10 mM Tris-HCl pH 7.5, 1 mM EDTA),  
2 X W/B buffer (1 M NaCl, 10 mM Tris-HCl pH 7.5, 1 mM EDTA),

ABgene: Reddy-Mix  
AnalaR: 100 % glycerol  
BDH: PCR-grade water  
Boehringer-  
Mannheim: *Sau3AI*  
DynaL: M2-80 Dynabeads: Streptavidin-coated magnetic beads  
FMC  
BioProducts: NuSieve GTG agarose  
Melford  
Laboratories: X-gal, IPTG  
Millipore: Microcon YM-100 / YM-50 spin column  
MWG-  
Biotech: Oligonucleotide A, PO<sub>4</sub>-Oligonucleotide B, 3'-biotinylated oligonucleotide,  
50 X UV-safe TAE buffer  
New England  
Biolabs: *AhaI*, *HaeIII*  
Promega: *RsaI*, 100 bp PCR ladder, 2 X T4 DNA Rapid ligation buffer, T4 DNA ligase,  
pGEM<sup>®</sup>-T vector, JM109 High Efficiency Competent Cells  
Qiagen: QIAquick gel extraction kit, MinElute PCR purification kit  
Sigma: ethidium bromide, ampicillin, N,N'-dimethyl-formamide

**Suggested Timetable:**

		<b>To do:</b>	<b>Approx. time:</b>	<b>To make:</b>
<b>Day 1</b>	<b>am</b>	Dilute oligonucleotides, Quantify genomic DNA.	30 min 45 min	2 % agarose gel, 20 X SSC, 6 X SSC,
	<b>pm</b>	2. Digest genomic DNA, Check digestion on 2 % agarose gel.	3 hr 45 min	2 X SSC, 1 X SSC, W/B Buffer.
<b>Day 2</b>	<b>am</b>	1. Anneal oligonucleotides,	1.5 hr	
	<b>pm</b>	2. Ligate adaptors to digested DNA.	30 min (o/night)	1.8 % NuSeive gel
<b>Day 3</b>	<b>am</b>	3. Size-select 400-1,000 bp fragments on 1.8 % Nuseive gel, Purify fragments from NuSeive gel,	1.5 hr 45 min	2 % agarose gel, Ampicillin.
	<b>pm</b>	PCR-test ligation success and size, Check PCR on 2 % agarose gel.	30 min (3 hr) 45 min	Agar plates, SOC medium.
	<b>Day 4</b>	<b>am</b>	4. Enrich and capture DNA.	2 hr
<b>Day 5</b>	<b>am</b>	4xxi. PCR enriched DNA fragments.	30 min (3 hr)	IPTG.
	<b>pm</b>	5. Clean-up PCR, Test PCR success on 2 % agarose gel,	45 min 45 min	X-gal.
		6a. Ligation into pGEM <sup>®</sup> -T vector.	30 min (o/night)	
<b>Day 6</b>	<b>am</b>	7b. Transformation into <i>E.coli</i> , Incubate competent cells.	1 hr 3 hr	IPTG/X-gal plates.
	<b>pm</b>	Plate out colonies and incubate overnight.	45 min	96-well 2 % agarose gel.
<b>Day 7</b>	<b>am</b>	8a. Pick positive clones, (9. PCR-screen positive clones).	upto 8 hr	30 % glycerol/LB, (Stabs).
	<b>pm</b>	(plate-out microsatellite-containing clones)	1 hr	
<b>Day 8</b>	<b>am</b>	Construct back-up/storage library, (Prepare stabs).	upto 5 hr 1 hr	
	<b>pm</b>			

**Notes:**

- Steps 4xxi-6a must be completed during the same day.

**1. Adaptor preparation:**

(i) Mix the following components in a 0.5 ml microcentrifuge tube:

10  $\mu$ l            **Oligonucleotide A** (200 pmol/ $\mu$ l)

10  $\mu$ l            **Oligonucleotide B** (200 pmol/ $\mu$ l)

(ii) Mix oligos by gentle pipetting or vortexing, then spin briefly.

(iii) Denature oligo mixture by heating at 80°C for 2-5 min on a thermocycler.

(iv) Allow oligos to anneal by leaving to cool at room temperature for 1 hour.

(v) Add 60  $\mu$ l of **PCR-grade water** (BDH).

[i.e. make up to a total of 80  $\mu$ l]

(vi) Store adapter mixture at -20°C.

Sticky end adaptor (*Sau3A1*-specific).

	5'		3'	
<b>Oligo A</b>		GGC CAG AGA CCC CAA GCT TCG		[21-mer]
<b>Oligo B</b>	PO <sub>4</sub> -	GAT CCG AAG CTT GGG GTC TCT GGC C		[25-mer]
Adaptor		GGC CAG AGA CCC CAA GCT TCG		
		CCG GTC TCT GGG GTT CGA AGC CTA G - PO <sub>4</sub>		[ <i>Sau3A1</i> restriction site]

Blunt end adaptor (e.g. *HaeIII* / *AluI* / *RsaI*).

	5'		3'	
<b>Oligo A</b>		CTC TTG CTT ACG CGT GGA CTA		[21-mer]
<b>Oligo B</b>	PO <sub>4</sub> -	TAG TCC ACG CGT AAG CAA GAG CAC A		[25-mer]
Adaptor		CTC TTG CTT ACG CGT GGA CTA		
		A CAC GAG AAC GAA TGC GCA CCT GAT - PO <sub>4</sub>		

**Notes:**

- Dilute oligonucleotides with PCR-grade water, leave to dissolve at 4°C for 2 hours to overnight, then make aliquots. Store oligonucleotide aliquots and stocks at -20°C.
- The adaptor consists of two oligonucleotides with complementary sequences that are heated and allowed to anneal together.
- Oligo B has a phosphate-group (PO<sub>4</sub>) at the 5' end and an extra four bases at its 3' end.

## 2. *Digestion of genomic DNA and adapter ligation:*

The importance of using high-quality genomic DNA for the construction of a genomic library *cannot* be emphasised enough - problems with sequencing likely happen if poor-quality DNA is initially used. It is also recommended that a test digest with several different enzymes is performed on *ca.* 1 µg of DNA to determine the restriction enzyme that provides the best size-distribution of DNA fragments.

All frozen components should be thawed on ice. Always add the restriction enzyme or T4 DNA ligase *last*, removing it from -20°C when needed and return the enzyme to -20°C *immediately*.

(i) Mix the following components in a 0.5 ml microcentrifuge tube:

. µl	genomic DNA	[9 µg]
9 µl	10 X enzyme buffer	[final concentration 1 X]
. µl	restriction enzyme	[40 U]
. µl	PCR-grade water (BDH)	
<b>90 µl</b>	<b>total volume</b>	

(ii) Incubate at 37°C for 1.5 - 2 hr.

(iii) Inactivate the enzyme by heating it to 65°C for 5-10 min.

(iv) To the digest add:

100 µl	2 X T4 DNA ligase buffer	[final concentration 1 X]
. µl	PCR-grade water (BDH)	
2.0 µl	adaptor	[25 pmol/µl]
2.0 µl	T4 DNA ligase (Promega)	[40 U]
<b>200 µl</b>	<b>total volume</b>	

(v) Mix gently by pipetting and incubate at either (a) 4°C overnight or (b) room temperature for 1 hr.

(vi) Inactivate the enzyme by heating to 65°C for 5-10 min.

(vii) If several samples are being digested (*e.g.* with more than one restriction enzyme) then the digested DNA can be pooled and concentrated using a **Microcon YM-50 spin column** (Millipore) if the total volume is too great.


(viii) Test the success of the digestion on a 2 % agarose gel made with **ethidium bromide** (Sigma) at a final concentration of 0.5 µg/ml.

### Notes:


- Suitable restriction enzymes include:
  1. *Sau3A1* for the sticky-end adaptor.
  2. *HaeIII*, *AluI* and *RsaI* for blunt-end adaptor.
- All restriction enzymes used in this protocol recognise a 4 bp restriction site - on average they are expected to cleave DNA every 256 bp.
- Using several blunt-end restriction enzymes and pooling the digestion product may increase the number of different fragments that contain microsatellites.

### 3. *Size-selection and PCR of adaptor-ligated DNA:*

- (i) Adaptor-ligated DNA is run (*ca.* 1 hour at 60 v) alongside a **100 bp PCR ladder** (Promega) in a 1.8 % **NuSieve GTG** (FMC BioProducts) agarose gel prepared with **UV-safe TAE buffer** (MWG Biotech) and **ethidium bromide** (Sigma) at a final concentration of 0.5 µg/ml.
- (ii) Visualise gel using a uv-transilluminator (**Figure 1**) and quickly (< 20 s) excise the 400-1,000 bp fragments (**Figure 2**) using the blunt end of a sterile scalpel blade (or a microscope cover-slip). Put the gel slices into a *pre-weighed* 1.5 ml centrifuge tube.



**Figure 1.** Digested genomic DNA alongside DNA ladder



**Figure 2.** Removal of 1,000-400 bp fragments.

- (iii) Purify the DNA fragments using a **QIAquick gel extraction kit** (Qiagen):
- (iv) Add 3 volumes of buffer QG to 1 volume of gel (100 mg ~ 100 µl).  
[maximum of 400 mg per QIAquick column]
- (v) Incubate at 50°C for 10 min (or until the gel slice has completely dissolved). *Briefly* (1-2 sec) vortex the tube every 2-3 min to help dissolve the sample.
- (vi) After the gel has dissolved check that the mixture is yellow.  
[If the colour is orange/violet then add 10 µl of 3 M sodium acetate pH 5.0 and mix - the mixture will turn yellow]
- (vii) Place QIAquick spin column in a 2 ml collection tube, apply sample to QIAquick column and centrifuge at room temperature for 1 min.
- (viii) Discard flow-through and replace QIAquick column in the same 2 ml collection tube.
- (ix) [optional] Add 0.5 ml buffer QG to QIAquick column and centrifuge for 1 min.
- (x) (Discard flow through) Add 0.75 ml buffer PE to QIAquick column, let stand for 2-5 min and centrifuge for 1 min.
- (xi) Discard flow-through and centrifuge the QIAquick for an additional 1 min at 13,000 rpm.

#### *Notes:*

- A 1.8 % agarose gel contains 0.54 g agarose in 30 ml buffer.
- Use a large-well comb so that enough digested DNA can be loaded into a single well.
- To prevent damage to the (*expensive*) uv-transilluminator keep the gel in the (*cheap*) casting tray when excising the fragments.
- 400-1,000 bp fragments are excised since fragments over *ca.* 1,000 bp long are unlikely to be completely sequenced, and sequences smaller than 400 bp are rarely long enough to provide good primers.

### 3. *Size-selection and PCR of adaptor-ligated DNA (cont.):*

(xiv) Place the QIAquick column into a clean 1.5 ml centrifuge tube.

(xv) Elute the DNA fragments in 30  $\mu$ l of 10 mM Tris-HCl (pH 8.5): apply elution buffer to centre of QIAquick column, let stand for 1 min and centrifuge at 13,000 rpm for 1 min.

[Buffer EB in Qiagen kits]

(xvi) Check the yield of recovered DNA on a 2 % agarose gel.

[approximately 4  $\mu$ g of DNA is expected]

(xvii) Pool all recovered DNA samples and reduce volume to *ca.* 10  $\mu$ l using a Microcon YM-100 (or YM-50) spin column (Millipore).

(xviii) Perform a PCR to test the success of the ligation of adaptors to the digested DNA; combine the following:

1.25 $\mu$ l	Oligonucleotide A (200 pmol/ $\mu$ l)
10.25 $\mu$ l	PCR-grade water (BDH)
1 $\mu$ l	size-selected DNA
12.5 $\mu$ l	Reddy-Mix (ABgene)
<b>25 <math>\mu</math>l</b>	<b>total volume</b>

(xix) Cycle the PCR mix on the following program:

Step 1	95°C	5 min
Step 2	95°C	50 s
Step 3	56°C	1 min
Step 4	72°C	2 min
Step 5	go to step 2 thirty times	
Step 6	72°C	10 min
Step 7	6°C	hold

[approximate time 2 hours]

(xx) Check PCR by running 5  $\mu$ l on a 2 % agarose gel alongside a 100 bp DNA ladder (Promega). Successful ligation of adaptors to digested DNA is indicated by a smear between 400-1,000 bp.

#### Notes:

- A 2 % agarose gel contains 0.6 g agarose in 30 ml buffer.
- Ensure that the elution buffer is dispensed directly onto the QIAquick membrane for complete elution of bound DNA.
- Elution efficiency is dependent on pH: maximum efficiency is achieved between pH 7.0 and 8.5.
- New Qiagen columns are available that elute in a 10  $\mu$ l volume.

#### 4. Capture of microsatellite-containing DNA fragments:

- (i) In a 1.5 ml microcentrifuge tube, wash 100  $\mu$ l of streptavidin-coated magnetic beads (10 mg/ml) (M-280 Dynabeads, Dynal) in 100  $\mu$ l of 1 X Washing/Binding (W/B) buffer (1 M NaCl, 10 mM Tris-HCl pH 7.5, 1 mM EDTA).
- (ii) Repeat (i) above.
- (iii) Re-suspend the beads in 200  $\mu$ l of 2 X W/B buffer.
- (iv) Add 2  $\mu$ l (200 pmol) of 3'-biotinylated oligonucleotide (at 100 mM).
- (v) Make the total volume to 400  $\mu$ l with PCR-grade water (BDH).
- (vi) Incubate the sample at room temperature for 30 min, with gentle agitation every 10 min.
- (vii) Wash once in 400  $\mu$ l of 1 X W/B and twice in 400  $\mu$ l of 6 X SSC.
- (viii) Re-suspend the beads in 50  $\mu$ l 6 X SSC and incubate at probe-specific hybridisation temperature (Table 1).
- (ix) In a separate tube add:

10 $\mu$ l	digested/adaptor-ligated/size-selected DNA	
. $\mu$ l	Oligonucleotide A (20 pmol)	
15 $\mu$ l	20 X SSC	[final concentration 6 X]
. $\mu$ l	PCR-grade water (BDH)	
<b>50 <math>\mu</math>l</b>	<b>total volume</b>	

- (x) Mix by gentle pipetting and briefly spin.
- (xi) Set up the following program on a thermal cycler:

Step 1	95°C	5 min	
Step 2	60°C	30 min	[probe-specific hybridisation temperature]
Step 3	70°C	2 hr	

- (xii) Incubate the DNA mixture (step ix) at 95°C for 10 min.
- (xiii) As the temperature ramps down to 60°C add the contents of the re-suspended bead mixture (steps i-viii) to the DNA sample, mix by gentle pipetting and incubate at 60°C for 30 min with gentle agitation every 5 min.
- (xiv) Remove beads using a magnet and quickly remove the supernatant. Re-suspend the beads in 100  $\mu$ l of 2 X SSC and transfer to a 1.5 ml microcentrifuge tube.  
[It is essential that the beads do not dry out between washes]
- (xv) Wash the bead mixture four times in 1 ml of 2 X SSC, incubating the sample for 5 min at room temperature before removing the 2 X SSC.
- (xvi) Wash the bead mixture four times in 1 ml of 1 X SSC, incubating the sample for 5 min at room temperature before removing the 1 X SSC.

#### Notes:

- To wash the Dynabeads<sup>®</sup> (i.e. remove supernatant) place the tube next to a magnetic stand (or magnetic stirrer) for 2-3 min to allow the beads to migrate to the side of the tube. Using a pipet, remove the supernatant without moving the tube or touching the beads. Next remove the tube from the magnet and add buffer as required; the beads may need resuspension by gentle pipetting. Repeat the steps as required.
- Removing beads is achieved by placing the tube next to a magnet for 2-3 min so that the beads are separated from the supernatant.



#### 4. *Capture of microsatellite-containing DNA fragments (cont.):*

- (xvii) Re-suspend the bead mixture in 100  $\mu$ l of **1 X SSC** and aliquot into four 25  $\mu$ l samples. Add 250  $\mu$ l **1 X SSC** to each aliquot and incubate at the probe-specific temperature for 10 min.  
 (xviii) Quickly remove the supernatant and rinse for 30 s at room temperature in 400  $\mu$ l **1 X TE**.  
 (xix) Rinse for 30 s at room temperature in 400  $\mu$ l **50 mM NaCl**.  
 (xx) Re-suspend each aliquot in 50  $\mu$ l **PCR-grade water** (BDH).

[final bead concentration is 5  $\mu$ g/ $\mu$ l]

- (xxi) Set up the following PCR to increase the quantity of enriched DNA:

25 $\mu$ l	<b>2 X Reddy-Mix</b> (ABgene)
3 $\mu$ l	<b>Oligonucleotide A</b> (10 pmol/ $\mu$ l)
8 $\mu$ l	bead suspension (40 $\mu$ g)
14 $\mu$ l	<b>PCR-grade water</b> (BDH)
<b>50 <math>\mu</math>l</b>	<b>total volume</b>

- (xxii) Run the PCR on the following program:

Step 1	95°C	3 min	
Step 2	95°C	30 s	
Step 3	60°C	30 s	[probe-specific annealing temperature]
Step 4	72°C	45 s	
Step 5	go to step 2 five times		
Step 6	92°C	30 s	
Step 7	60°C	30 s	[probe-specific annealing temperature]
Step 8	72°C	55 s	
Step 9	go to step 6 thirty times		
Step 10	72°C	30 min	
Step 11	4°C	hold	[approximate time 3 hours]

- (xxiii) Run 2-5  $\mu$ l of PCR product alongside a **PCR marker** (Promega) on a **2 % agarose gel** at 100 V for 15-20 min. A smear between 0.4 - 1 Kbp indicates a successful capture and PCR.

**Table 1.** Hybridisation, washing and PCR temperatures for some oligonucleotides.

Probe	T <sub>m</sub> / °C	Hyb / °C	PCR T <sub>a</sub> / °C	Probe	T <sub>m</sub> / °C	Hyb / °C	PCR T <sub>a</sub> / °C
(GT) <sub>15</sub>	82	75	65	(AAT) <sub>8</sub>	56	50	40
(GT) <sub>12</sub>	77	70	60	(CAG) <sub>8</sub>	84	78	68
(AAAG) <sub>6</sub>	77	67	57	(GAA) <sub>8</sub>	80	70	60
(GTG) <sub>8</sub>	94	84	64	(AAC) <sub>8</sub>	70	64	54

#### *Notes:*

- Annealing temperature of probes may be calculated using the oligo calculator at: <http://www.basic.nwu.edu/biotools/oligocalc.html>.
- For PCR use an annealing temperature 5-10°C below the annealing temperature calculated above.
- Proceed to the PCR clean-up and ligation steps immediately (or the next day at the latest!).

**5. PCR clean up:**

(i) Add 5 volumes of buffer PB to 1 volume of the PCR reaction and mix.

[It is not necessary to remove mineral oil]

(ii) Place a MinElute column in a 2 ml collection tube.

(iii) Apply the sample to the MinElute column and centrifuge for 1 min at 13,000 rpm.

(iv) Discard flow-through. Place the MinElute column back in the same 2 ml collection tube.

(v) Add 750  $\mu$ l buffer PE to the MinElute column and centrifuge for 1 min at 13,000 rpm.

(vi) Discard flow-through and place the MinElute column back in the same collection tube.

Centrifuge the column for 1 min at 13,000 rpm.

[Residual ethanol from buffer PE will not be removed unless the flow-through is discarded before this step]

(vii) Place the MinElute column in a clean 1.5 ml centrifuge tube.

(viii) To elute the DNA, add 10  $\mu$ l buffer EB (or 10 mM Tris-HCl, pH 8.5) to the centre of the membrane, let the column stand for 1 min and centrifuge for 1 min at 13,000 rpm.

---

**Notes:**

**6a. Ligation of enriched DNA into pGEM<sup>®</sup>-T vector.**

- (i) Briefly centrifuge the **pGEM<sup>®</sup>-T vector** to collect the contents at the bottom of the tube.  
(ii) Vortex the **2 X Rapid ligation buffer** vigorously before each use.  
(iii) Using 0.5 ml centrifuge tubes with low DNA-binding capacity set up the following ligation reaction:

5 µl	<b>2 X Rapid Ligation Buffer</b> (Promega)
1 µl	<b>pGEM<sup>®</sup>-T vector</b> (50 ng) (Promega)
x µl	PCR product ( <i>see section 6b below</i> )
1 µl	<b>T4 DNA ligase</b> (3 Weiss units/µl)
x µl	<b>PCR-grade water</b> (BDH)
<b>10 µl</b>	<b>total volume</b>

- (iv) Mix the reactions by pipetting.  
(vi) Incubate the reaction for either  
(a) 1 hr at room temperature, or  
(b) overnight at 4°C. [this provide the maximum number of transformants]

**6b. Optimising insert:vector molar ratios.**

- (i) **pGEM<sup>®</sup>-T vector** was optimised using a 1:1 molar ratio of insert DNA to vector - ratios of 1:8 to 8:1 can be used.  
(ii) Ratios of 3:1 to 1:3 are typically used initially.  
(iii) **pGEM<sup>®</sup>-T vector** is approximately 3 Kb long.  
(iv) To calculate the appropriate amount of PCR product (insert) to use in the ligation reaction use the following equation:

$$\frac{\text{ng vector} \times \text{Kb size of insert}}{\text{Kb size of vector}} \times \text{insert:vector molar ratio} = \text{ng of insert}$$

**Notes:**

- Use the Promega T4 DNA ligase supplied with the pGEM<sup>®</sup>-T vector kit.
- 2X rapid ligation buffer contains ATP which degrades during temperature fluctuations; avoid multiple freeze-thaw cycles by making single-use aliquots.
- **IMPORTANT:** Vortex the 2X rapid ligation buffer before each use.
- Longer incubation times will increase the number of transformants.
- For overnight ligations, float the tubes in polystyrene on a bath of room-temperature water and place the bath in 4°C refrigerator.

### 7a. Preparation of agar plates.

- (i) **LB/ampicillin plates**: Melt 400 ml **LA** in microwave (at *defrost* for 15-20 min). Allow **LA** to cool to 50°C before adding **ampicillin** to final concentration of 100 µg/ml. Pour medium (*ca.* 30-35 ml) into petri-dishes and allow to harden. Store plates at 4°C (for up to 1 month).
- (ii) **LB/ampicillin/IPTG/X-gal plates**: Spread 100 µl of 0.1 M **IPTG** and 20 µl of 50 mg/ml **X-gal** over the surface of an **LB/ampicillin** agar plate; allow to absorb for at least 30 min at 37°C prior to use.
- (ii) **LB/ampicillin/IPTG/X-gal plates**: [*Optional*] Make **LB/ampicillin plates** as above but add 0.5 mM **IPTG** and 80 µg/ml **X-gal** before pouring plates.

### 7b. Transformation of pGEM<sup>®</sup>-T vector into *E. coli*.

- (i) Prepare **LB/ampicillin/IPTG/X-gal** plates as described above. Equilibrate plates to room temperature prior to plating.
- (ii) Centrifuge tubes containing ligation reaction to collect contents at bottom of tube.
- (iii) Add 2 µl of each ligation reaction to a 1.5 ml centrifuge tube on ice.
- (iv) Remove tube(s) of frozen **JM109 High Efficiency Competent Cells** (from -70°C storage) and place in an ice bath until *just* thawed (*ca.* 5 min). Mix the cells by *gently* flicking the tube.
- (v) *Carefully* transfer 50 µl of cells into each tube prepared in step 2.
- (vi) *Gently* flick the tubes to mix and place them on ice for 20 min.
- (vii) Heat-shock the cells for 45-50 sec in a water bath at *exactly* 42°C - **DO NOT SHAKE**.
- (viii) Immediately return the tubes to ice for 2 min.
- (ix) Add 950 µl room temperature **SOC medium** to the tubes containing cells transformed with ligation reactions. [LB broth may be substituted, but colony number may be lower]
- (x) Incubate for 1.5 hr at 37°C with shaking at approximately 150 rpm.
- (xi) Plate 100 µl of each transformation culture onto **LB/ampicillin/IPTG/X-gal** plates.
- (xii) Incubate the plates overnight (16-24 hr) at 37°C.
- (xiii) [*optional*] incubation of plates at 4°C (after 37°C incubation) helps develop the blue/white colour of colonies.

#### Notes:

- Dissolve ampicillin in sterile distilled water to a final concentration of 50 mg/ml; store as aliquots at -20°C. Adding ampicillin to **LA** before it cools below 50°C destroys the ampicillin.
- Colonies generally grow better on plates that have been prepared a few days beforehand.
- Be *extremely* careful with competent cells - any vigorous treatment will destroy them.
- The use of 17x100 mm polypropylene tubes (Falcon # 2059) usually increases transformation efficiency - tubes from some manufacturers bind DNA and decrease colony number.

**8a. Working library construction:**

- (i) Add **ampicillin** to **LB** to a final concentration of 100 µg/ml.
- (ii) Pipette 100 µl of **LB/ampicillin** into each well of as many 96-well plates as required.
- (iii) Using a sterile **toothpick** or a **200 µl pipet tip** (one for each colony), pick a colony containing an insert and swirl tip into a well; repeat until all colonies have been picked making sure that only one colony is added to each well.
- (iv) Incubate 96-well plates at 37°C for 3-4 hours.
- (v) Store plates at 4°C for up to 3 weeks.

**8b. Long-term storage library construction:**

- (i) Prepare glycerol stock medium by adding 30 ml **100 % glycerol** (AnalaR) to 70 ml **LB**.  
[final glycerol concentration of 30 %]
- (ii) Autoclave 30 % glycerol stock.
- (iii) Add **ampicillin** to **LB** to a final concentration of 100 µg/ml.
- (iv) Pipette 100 µl of **LB/ampicillin** into each well of as many 96-well plates as required.
- (iv) Using sterile 1-10 µl pipet tips transfer 10 µl of the working library into corresponding wells of the back-up plates.
- (v) Incubate 96-well plates at 37°C for 3-4 hours.
- (vi) Add 100 µl of 30 % glycerol stock to each well.
- (vii) Store at -80°C. [Bacterial cultures may be stored indefinitely in glycerol cultures]

**8c. Stab construction:**

- (i) Melt **LA** in microwave (15-20 min on defrost).
- (ii) Add **ampicillin** to **LB** to a final concentration of 100 µg/ml when **LB** is at 50°C or lower.
- (iii) Use 2 ml plastic vials with screw caps fitted with rubber gaskets.
- (iv) Add 1.5 ml molten **LA/ampicillin** and let cool with vial caps loose.
- (v) Tighten vial caps when the agar has set.
- (vi) To store bacteria, pick a single, well isolated colony with a sterile loop.
- (vii) Stab the loop through the agar to the bottom of the vial.
- (viii) Replace and tighten the cap. Label the vial and cap.
- (ix) Store the vial in the dark at room temperature.

---

**Notes:**

- White colonies generally contain inserts, but see notes in the Promega pGEM<sup>®</sup>-T vector cloning kit manual.

**9a. PCR screening for microsatellite-containing clones.**

(i) Prepare the following master mix:

5 $\mu$ l	Oligo A (10 pmol/ $\mu$ l)
2.5 $\mu$ l	(GT) <sub>15</sub> oligonucleotide (10 pmol/ $\mu$ l)
2.5 $\mu$ l	PCR-grade water (BDH)
10 $\mu$ l	2 X Reddy Mix (ABgene)
<b>20 <math>\mu</math>l</b>	<b>total volume</b>

(ii) Pipette 20  $\mu$ l of the PCR master mix into each well of 96-well PCR plate.


(iii) At the same time as **stage 8a** place the toothpick/pipet tip into a well corresponding to the library plate and swirl to release some cells.

(iii) Run the PCR on the following PCR program (SCRGT):

Step 1	95°C	3 min	
Step 2	95°C	30 s	
Step 3	60°C	30 s	[probe-specific annealing temperature]
Step 4	72°C	45 s	
Step 5	go to step 2 five times		
Step 6	92°C	30 s	
Step 7	60°C	30 s	[probe-specific annealing temperature]
Step 8	72°C	55 s	
Step 9	go to step 6 thirty times		
Step 10	72°C	30 min	
Step 11	4°C	hold	[approximate time 3 hours]

(iv) Load 5-10  $\mu$ l of PCR product alongside a **PCR marker** (Promega) on a **2 % agarose gel** containing **ethidium bromide** (Sigma) at a final concentration of 0.5  $\mu$ g/ml. Run the gel at approximately 100 V for about 15-20 min.

(v) A double banded PCR product indicates that the vector contains a microsatellite-containing insert (**Figure 3**).



**Figure 3.** PCR screening of clones. Blue circles indicate a microsatellite.

Notes:

**APPENDIX: Commonly used buffers and solutions.**

- Ampicillin:** Dissolve 50 mg in 1 ml H<sub>2</sub>O and filter sterilize. Store in 1 ml aliquots at -20°C.
- 0.5 M EDTA:** Add 186.1 g disodium ethylenediaminetetra-acetate.2 H<sub>2</sub>O (EDTA) to 800 ml H<sub>2</sub>O. Stir vigorously on a magnetic stirrer. Adjust pH to 8.0 with NaOH (*ca.* 20 g NaOH pellets). Sterilise by autoclaving.
- Ethidium Bromide:** Add 1 g ethidium bromide to 100 ml of H<sub>2</sub>O and stir on magnetic stirrer for several hours. wrap the container in aluminium foil and stor at room temperature. The solution concentration is 10 mg/ml.
- 0.1 M IPTG:** Dissolve 1.2 g IPTG in water to 50 ml final volume. Filter sterilise and store at -20°C.
- LA:** Prepare LB as below. Just before autoclaving add 15 g (per litre) of bacto-agar.
- LB medium:** To 950 ml H<sub>2</sub>O add 10 g bacto-tryptone, 5 g bacto-yeast extract, 10 g NaCl. shake the solution until solutes dissolve and adjust the volume to 1 litre. Sterilise by autoclaving.
- 2 M Mg<sup>2+</sup>:** Mix 20.33 g MgCl<sub>2</sub>.6H<sub>2</sub>O and 24.65 g MgSO<sub>4</sub>.7H<sub>2</sub>O with distilled water to 100 ml; filter sterilise and store at 4°C.
- 1 M NaCl:** Dissolve 58.44 g of NaCl in 1 litre of H<sub>2</sub>O.
- 50 mM NaCl:** Dissolve 2.922 g of NaCl in 1 litre of H<sub>2</sub>O or dilute 1 M NaCl.
- SOC: :** Dissolve 2.0 g Bacto<sup>®</sup>-tryptone, 0.5 g Bacto<sup>®</sup>-yeast extract, 1 ml 1 M NaCl, 0.25 ml 1 M KCl, 1 ml 2 M Mg<sup>2+</sup> stock, 1 ml 2 M filter-sterilised glucose. Store at 4°C.
- 20 X SSC:** Dissolve 175.3 g NaCl and 88.2 g sodium citrate in 800 ml of H<sub>2</sub>O. Adjust to pH 7.0 with a few drops of a 10 N solution of NaOH. Make volume to 1 litre with H<sub>2</sub>O and autoclave.
- 50 X TAE:** Dissolve 242 g Tris base, 57.1 ml glacial acetic acid and 100 ml 0.5 M EDTA (pH 8.0) in 1 litre of water. *Note:* this solution is 50 X - dilute 49:1 for working concentration.
- 1 M Tris:** Dissolve 121.1 g of Tris base in 800 ml of H<sub>2</sub>O. Adjust pH to desired value (*e.g.* 8.5 or 7.4) with HCl at room temperature; adjust the final volume to 1 litre with H<sub>2</sub>O. Aliquot and sterilise by autoclaving. Other Tris buffers can be prepared by adjusting an aliquot of 1 M Tris with an appropriate volume of water.
- 1 X W/B buffer:** Dissolve 5.844 g NaCl, 37 mg EDTA in 50 ml 20 mM Tris-HCl pH 7.5. Make volume up to 100 ml with distilled water. The final concentration of the solution is: (1 M NaCl, 10 mM Tris-HCl pH 7.5, 1 mM EDTA).
- 2 X W/B buffer:** Dissolve 11.688 g NaCl, 37 mg EDTA in 50 ml 20 mM Tris-HCl pH 7.5. Make volume up to 100 ml with distilled water. The final concentration of the solution is: (2 M NaCl, 10 mM Tris-HCl pH 7.5, 1 mM EDTA).
- X-gal:** Dissolve 100 mg 5-bromo-4-chloro-3-indolyl-β-D-galactosidase in 2 ml N,N'-dimethyl-formamide. Cover with aluminium foil and store at -20°C. The final concentration of X-gal is 50 mg/ml.